## PTD-100 User manual

## / Components /


Product



Additional Cover

## / Product Contents /

## / How to Use /

## / Installing the application /

Install Thermo Capsule app on your smartphone

Please search for thermo capsule in Google Play or Apple APP Store and install it.

Running application

Run the application by touching the icon of thermo capsule (thermo capsule application runs on Android 4.2 or higher or on IOS 7.0 or higher).

## / Connecting Thermo Capsule /

Connect the capsule to the earphone jack of your smartphone.

Connect the capsule to the earphone jack of your smartphone.

Connect the capsule as shown on the right and open the slide cap. (When the connection is made normally, the icon at the bottom left corner of the screen lights on green).

If the earphone jack is located at the bottom of the smartphone, more convenient measurement is possible by touching the icon (🔄) at the bottom center of the screen.


## / Body Temperature /

The body temperature can be measured immediately.

When the click button on the APP screen is touched, after having connected the body temperature capsule, the body temperature is measured.


If the click button is pressed longer, the user's own body temperature is measured conveniently by stating the countdown. Check the body temperature displayed at the top and save it by touching the Save button at the bottom right corner if you want to save it.

## / Statistics /

The measured temperature.

If thermo capsule APP is used, body temperature can be measured and saved. The saved body temperatures can be viewed on a graph.


## / Object Temperature /

Not only body temperature but also object's temperature can be measured accurately.

Specific object's emissivity can be measured more accurately if this specific object is selected before the measurement process. The measurement of emissivity can be applied by selecting the icon at the bottom of the screen.

Choose the EMS mode for more diverse objects.


The air temperature should be measured after having closed the slide cap of the capsule. For accurate measurement, expose the capsule for more than 10 min in the air.

## / Alarm /

Alarms can be set in order to remind the user to measure the temperature regularly.


The alarm works like a timer and it sounds after the chosen time has elapsed.

## / Search /

The Facilities such as neighbouring hospitals, pharmacies or emergency rooms, etc can be immediately found.

You can search and find nearby facilities, for which the device uses Google maps (it may vary depending on location). In case of very high body temperature, the condition may be considered as a medical emergency, for which visit at nearby medical facility to receive medical diagnosis is recommended.

Product Specification	
Product Name	THERMO CAPSULE
Name of Item	Skin infrared thermometer
Model	PTD-100
Item Authorization No.	JEIN 15-4136
Purpose	Instrument for measuring body temperature (and others) by using infrared thermography of skin (like forehead, etc.)
Measurement range	Body temperature : 32~42°C Accuracy : ±0.2°C (36~39°C) / other accuracy : ±0.3°C
Measurement time	Within 5 seconds
Measuring method	Infrared radiation energy measurement
Rated power	DC 3.0V (1.5DC, 250mAh x 2)
Weight	13.7g (main body 9.2g, Cover 3.8g)
Communication method	Ø3.5 earphone jack (quadripolar)/ CTIA
Component	1 Set + 1 additional cover
Compatibility	Android 4.2 or above / iOS 7.0 or above

## / Precaution /

1. Dropping or shaking the product may damage the product.
2. Since the sensor area is very sensitive, please maintain it clean and not to be damaged for accurate measurement.
3. Do not expose the product to flammable materials and direct sunlight and store the product in a place without dust and pollution.
4. Keep it away from water.
5. It is recommended to store the product at room temperature. In case the product is not stored at room temperature, leave it for 10 minutes at room temperature before using it.
6. Clean the sensor areas carefully with a cotton swap, moistened in alcohol.
7. When measuring, please maintain 3~5cm of distance between the capsule and the object.
8. If the posture of the product is not right when measuring the body temperature, the measured value may be inaccurate.
9. Please measure the body temperature once the product is stabilized at the room temperature. After working out, showering, or being exposed to rain, the body temperature measurement may be inaccurate.
10. For accurate measurement, please measure the body temperature after having removed hair, cosmetics, etc. from the area to be measured.
11. For accurate measurement measure the temperature on both body lobes. The temperature values measured on the right and left side of the forehead may differ (from each other).
12. If the smartphone's battery is not sufficiently charged, the measurement may be inaccurate. Please use the product after charging the battery.
13. When measuring, place the sensor area close to the previously cleaned object and wait for the vibration of the measuring device, which announces that the measurement is completed. After this results can be seen on the smart phone screen.
14. If the person whose body temperature needs to be measured has just comes inside from the outside, or in case the measuring device is not at room temperature waiting at least 30 minutes before measuring is recommended.
15. Do not allow children to suck the product or its battery or touch the measuring instrument with wet hands.
16. Before the measurement, do not touch the slide cap or the sensor area. The temperature may not be measured accurately.
17. Repeat the measurement more than three times under following circumstances:
  - For 90 days old babies
  - For infants under 3 years only, the average temperature range is displayed
  - If the product is being used for the first time or the user has no previous experience of using any similar product
18. If you are sweating on your forehead, we recommend to measure the temperature behind your ear lobes. Please throw your hair back completely and hold the thermometer behind your ear lobe vertically. Please move, shortly the thermo-meter behind your ear lobes up and down.

## / Battery Replacement /

The battery can be replaced with the battery sold in market.

1. Remove the screw at the back of the main unit (red circle) with a precision screwdriver.
2. Separate the back of the main unit by pushing up the part as the red arrow shows on the picture above.
3. Replace the battery after disassembling front and back of the product. Assembly is done in reverse order.
4. The product contains two preinstalled batteries at 25mAh, LR41, which can be easily purchased.
5. Battery can be used about 1 year based on 15 times of measurement per day. (Body temperature measurement)

## \* Usage!

1. Can be used to measure the body temperature of babies.
2. Can be used to measure the air temperature.
3. Can be used to measure the temperature of different objects, for example, mobile phone, notebook, etc. If EMS mode is used, the temperature of the object can be measured more accurately.
4. Can be used to measure food and water temperature while cooking.
5. Can be used to measure bath temperature.
6. Can be used to measure the temperature of the milk for the baby.

## / Product Warranty /

Product Name	THERMO CAPSULE		
Purchase Date	Y	M	D
Purchasing Office			
Model Name	PTD-100		
Warranty Period	2 year warranty from the date of purchase		

**partron**

Web Manual Asia

www.partron.co.kr/manuals

**PARTRON**


22, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, Korea

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PARTRON Europe GmbH